

# Oscilloscopes for PDV: Getting More out of your Scope

Stephen Mueller


Applications Engineering Manager

[Stephen.Mueller@LeCroy.com](mailto:Stephen.Mueller@LeCroy.com)



**TELEDYNE LECROY**  
Everywhereyoulook™


# Oscilloscope Block Diagram



- ✦ Amplifiers and attenuators scale the input signals amplitude to match the input range of the Analog to Digital Converter (ADC)
- ✦ The ADC samples the analog signal and converts the voltage into a digital number

# Quantization Error

*4x attenuation  
of signal  
produces 4x  
larger  
quantization  
steps = 2 bits  
less resolution!*



Zooms of  
rising edges  
show  
quantization  
levels

# Enhanced Resolution

## Description


- ★ ERES increases vertical resolution by filter the waveform
- ★ Gain up to 3 bits of resolution in ½ bit steps

## Usage

- ★ Suitable for single-shot and repetitive waveforms
- ★ Reduces noise
- ★ Increases precision

## Disadvantages

- ★ Effective bandwidth decreased by 50% for every ½ bit gain



# 12-bit High Definition Technology – Teledyne LeCroy HDO6000



- ✦ Combination of
  - ✦ High Sample Rate 12-bit ADCs
  - ✦ High signal-to-noise input amplifiers
  - ✦ Low noise system architecture
- ✦ 16 times more resolution than traditional 8-bit scope
- ✦ Capture high frequency signals with 1GHz bandwidth
- ✦ Benefits
  - ✦ Clean, Crisp Waveforms
  - ✦ More Signal Details
  - ✦ Precise Waveform Measurements

# Introducing the Teledyne LeCroy HDO8000 Oscilloscope

## ★ *Further*

- ★ 8 analog input channels
  - ★ Ideal for high power and three-phase power electronics analysis
  - ★ Very useful for deeply embedded electronic/mechatronic systems

## ★ *Finer*

- ★ 12-bit HD4096 High Definition Technology
  - ★ “16x closer to perfect”

## ★ *Faster*


- ★ Up to 1 GHz
  - ★ Bandwidth for today's **and** tomorrow's technologies



*HD4096 12-bit technology is being deployed into 8 channels to meet the needs of fast growing applications*



# Traditional High Channel Count Setup: Coherent Modulation Measurements, 2009-2010



# LabMaster 10-100 Zi 100 GHz Real Time Bandwidth

- ✦ Starting at 25 GHz bandwidth
- ✦ The highest bandwidth: Up to 100 GHz
- ✦ The fastest sample rate: Up to 240 GS/s
- ✦ The highest channel count:
  - ✦ Up to 80 channels @ 36 GHz
  - ✦ Up to 40 channels @ 65 GHz
  - ✦ Up to 20 channels @ 100 GHz






# 60 Channel 10Zi System at NY Factory



# Track View


*Frequency is measured and plotted for each cycle*




*Frequency of Track plot reveals FM modulation rate of 100 kHz*

- ✦ Histograms reveal statistical distribution of values, but they contain no timing information
- ✦ Track View plots the value of a measurement versus time
- ✦ Track waveform values are calculated cycle by cycle
- ✦ Track waveform is time coincident with source waveform

# Track of Frequency Reveals Shape of Velocity Curve



# User-Defined Custom Measurements






***Custom parameters are selected from the measurement menu like any standard parameter. Supports MATLAB, C/C++, Visual Basic Script, VBA, MS Excel, and JavaScript***




# User-Defined Custom Math

**Custom math functions are selected from the math setup menu like any standard function. Script is entered in**




# In-Line Custom Math and Measurements








# MATLAB Spectrogram on Scope



# Create Custom Math Functions with Excel




Trace generated in  
Excel



Measurements on Excel-generated trace



# MATLAB Measurements



- ✦ Custom MATLAB parameter finds the half-life of the damped sine wave
- ✦ The value 3.149 is the number of cycles that have occurred when the signal reaches 50% of its peak amplitude